

## The EZee Fold portable power wheelchairs

The CH4075 EZee Fold is a high weight capacity folding portable electric wheelchair. With a 20" x 20" seat size and a 500 lb weight capacity it will fit larger users who fall within these specifications. Like all other EZee Fold wheelchairs, the CH4075 easily folds into a compact package that will fit in most vehicles. This is an excellent choice for people who want a portable folding electric wheelchair that can be used both indoors and out.

### CH4075 Bariatric Electric Folding Wheelchair

The CH4085 is the largest of our EZee Fold line of portable folding wheelchairs offering a 20" seat and weight capacity of 500 pounds it will take larger adults where they want to go while still providing a folding frame for easy storage and transport.

**Features:**

- 500 lb Capacity
- Right or Left Hand Drive
- Powered or Manual Modes
- 20 KM Range
- 8 KPH Maximum Speed

Length:	128 cm (50.4")
Width:	65.6 cm (25.8")
Height:	102.9 cm (40.5")
Folded Width:	44.8 cm (17.6")
Seat Width:	50.8 cm (20")
Seat Depth:	50.8 cm (20")
Seat Height:	54 cm (21.3")
Back Height:	45.7 cm (18")
Armrest Height:	Adjustable
Chair Weight:	38.2 kg (84 lb) w/battery
Weight Capacity:	227 kg (500 lb)
Caster Size:	20 cm (8")
Wheel Size:	30 cm (12")
Motors (2):	50 amp, brushless motor controller
Battery:	24 V, 20 Ah (3.2 kg/7 lb)
Charger Output:	5 Amps
Maximum Speed:	8 kph (5 mph)
Min. Turning Radius:	99 cm (39")
Cruising Range:	20 km (12.5 miles)
Net Weight:	25 kg (55 lb)
Gross Weight:	35 kg (77.2 lb)
Box Dimensions:	85 x 43 x 82 cm



Pictures and specifications are for reference only and may differ slightly.

