

The EZee Fold portable power wheelchairs

are the ideal solution for those who want to take a power wheelchair along with them on trips and excursions but don't want to have to rely on wheelchair vans and lifts. For added performance the EZee Fold wheelchairs can be fitted with 2 batteries to extend their range. Additional batteries must be ordered separately.

CH4050 EZee Fold Standard

Weighing approximately 51.4 lbs, including one battery, the EZee Fold Standard can easily be folded and lifted into most vehicles while still having a weight capacity of 264 lbs.

Features:

- Folds and opens in just 3 seconds.
- Strong and lightweight frame made of high grade aircraft aluminum alloy.
- Stands upright after folding for easy transport and storage.
- Li-ion batteries
- Flip up armrest for easier sliding transfers and to make tables more accessible
- Folding footrest for easier transfers
- With adjustable back belt and detachable back cushion
- With detachable seat cushion and washable seat cushion cover.

Seat Width:	460 mm (18.1")
Seat Height:	460 mm (18.1")
Seat Depth:	430 mm (17")
Backrest Height:	450 mm (17 3/4")
Net Weight:	23.1 kg (51.4 lbs.) with 1 battery
Weight Capacity:	120kgs (264 lbs)
Caster:	8" PU Tire
Rear Wheel:	10" PU Tires
Battery:	Lithium battery, 24 V 9 Ah
Max speed:	6.5 kph (4 mph)
Max. Climbing Grade:	10 Degrees
Charging period:	5 Hours per battery
Charger:	Universal charger
Endurance:	15 km (9.4 miles) / battery
Overall size:	870 x 560 x 910 mm (34.2" x 22" x 35.8")
Folded size:	560 x 300 x 660 mm (22" x 11.8" x 25.9")



CH4050 will hold either 1 or 2 batteries. All batteries are charged simultaneously on chair. Fixed push handle height.

Performance will be affected by user weight and battery condition.

CH4055 EZee Fold Swing Away (8" Wheels)

Weighing approximately 51.4 lbs, including one battery, the EZee Fold CH4055 can easily be folded and lifted into most vehicles while still having a weight capacity of 286 lbs.

Features:

- Folds and opens in just 3 seconds.
- Strong and lightweight frame made of high grade aircraft aluminum alloy.
- Stands upright after folding for easy transport and storage.
- Li-ion batteries
- Flip up armrest for easier sliding transfers and to make tables more accessible
- Folding footrest for easier transfers
- With adjustable back belt and detachable back cushion
- With detachable seat cushion and washable seat cushion cover.

Seat Width:	400 mm (16")
Seat Height:	460 mm (18.1")
Seat Depth:	430 mm (17")
Backrest Height:	450mm (17 3/4")
Net Weight:	23.1 kg (51.4 lbs.) with 1 battery
Weight Capacity:	120 kg (264 lbs.)
Caster:	8" PU Tire
Rear Wheel:	10" PU Tire
Battery:	Lithium battery, 24 V 9 AH
Max speed:	6.5 kph (4 mph)
Max. Climbing Grade:	10 Degrees
Charging period:	5 Hours per battery
Charger:	Universal charger
Endurance:	Up to 15 km (9.4 miles) /battery
Overall size:	960 x 589 x 935 mm (37.8" x 23.2" x 36.8")
Folded size:	560 x 300 x 660 mm (22" x 11.8" x 25.9")



Pictures and specifications are for reference only and may differ slightly.

The EZee Fold portable power wheelchairs

are the ideal solution for those who want to take a power wheelchair along with them on trips and excursions but don't want to have to rely on wheelchair vans and lifts. For added performance the EZee Fold wheelchairs can be fitted with 2 batteries to extend their range. Additional batteries must be ordered separately.

CH4060 EZee Fold HD 12" Rear Wheels

Weighing approximately 56 lbs, including one battery, the EZee Fold HD can easily be folded and lifted into most vehicles while still having a weight capacity of 350 lbs.

Features:

- Folds and opens in just 2 seconds.
- Strong and lightweight frame made of high grade aircraft aluminum alloy.
- Li-ion batteries.
- Flip up armrest for easier sliding transfers and to make tables more accessible.
- Folding footrest for easier transfers
- With adjustable back belt and detachable back cushion.
- With detachable seat cushion and washable seat cushion cover.

Seat Width:	457 mm (18")
Seat Height:	422 mm (16.6")
Seat Depth:	400mm (15 3/4")
Height adjustable backrest:	450 to 650 mm range (17.7" to 25.6")
Net Weight:	25.3 kg (56 lbs.) with 1 battery
Weight Capacity:	160 kg (350 lbs.)
Caster:	8" PU Tire
Rear Wheel:	12" Rubber/Pneumatic tire
Battery:	Lithium battery, 24V, 10AH
Max speed:	7.2 kph (4.5 mph)
Max. Climbing Grade:	12 Degrees
Charging period:	5 Hours per battery
Charger:	Universal charger
Endurance:	15 km (9.4 miles) battery
Overall size:	1000 x 640 x 890 mm (39.4" x 25.2" x 35")
Folded size:	640 x 335 x 740 mm (25.2" x 13.2" x 29.1")



CH4060 will hold either 1 or 2 batteries. All batteries are charged simultaneously on chair. Height adjustable push handle 450 mm to 650 mm (17.7" to 25.6"). Performance will be affected by user weight and battery condition.

CH4065 EZee Fold Swing Away (12" Wheels)

Weighing approximately 58.6 lbs, including one battery, the EZee Fold CH4065 can easily be folded and lifted into most vehicles while still having a weight capacity of 350 lbs.

Features:

- Folds and opens in just 2 seconds.
- Strong and lightweight frame made of high grade aircraft aluminum alloy.
- Li-ion batteries.
- Flip up armrest for easier sliding transfers and to make tables more accessible.
- Folding footrest for easier transfers
- With adjustable back belt and detachable back cushion.
- With detachable seat cushion and washable seat cushion cover.

Seat Width:	430 mm (17")
Seat Height:	507 mm (20")
Seat Depth:	430 mm (17")
Backrest Height:	441 mm (17.4")
Net Weight:	26.6 kg (58.6 lbs.)
Weight Capacity:	160 kg (350 lbs.)
Caster:	8" PU Tire
Rear Wheel:	12" Rubber Tire
Battery:	Lithium battery, 24V, 10 AH
Max speed:	7.2 kph (4.5 mph)
Max. Climbing Grade:	12 Degrees
Charging period:	5 hours per battery
Charger:	Universal charger
Endurance:	Up to 15 km (9.4 miles) / battery
Overall size:	954 x 629 x 948 mm (37.6" x 24.8" x 37.3")
Folded size:	629 x 346 x 752 mm (24.8" x 13.6" x 29.6")



Pictures and specifications are for reference only and may differ slightly.