

4-IN-1 MULTIFUNCTIONAL ELECTRIC LIFT TRANSFER, SHOWER, & COMMUNE CHAIR

CH3075

This versatile patient commode/shower transfer chair serves multiple functions: it can be used as a shower chair, wheelchair, commode chair, over-the-toilet commode, and transfer chair. The electric height adjustment feature facilitates effortless transfers to and from the bed, eliminating the need for caregivers to manually lift patients. This innovative design preserves the client's dignity compared to traditional sling transfer systems. Additionally, it can be used over the toilet for client toileting

SPECIFICATIONS

Material	Steel
Weight Capacity	135kg (300lbs)
Dimensions	131cm *54.5cm* 72.5cm
Chair weight	28.5kg (63lbs)
Seat height	16.5" - 29.5" (42cm-75cm)
Internal width	18" x 13.5" (45.5cm-48cm)
Seat width	18" x 13.5" (45.5cm-34cm)
Arm height	13" (34cm)
Back height	16" (42cm)
Padded arm length	13" (34cm)
Arm length	15" (38cm)

All measurements are approximate.



4-IN-1 MULTIFUNCTIONALITY TRANSFER ELECTRIC LIFT SHOWER/COMMUNE CHAIR


CH3075

BENEFITS

- Comfortable padded seat, back, & arms enhance user experience.
- Client transfers improve dignity by allowing clients to sit during transfers.
- No lifting required for PSWs or support workers, ensuring worker safety.
- Protective cover over the battery & electric lifting mechanism allows safe use in water environments.
- 180-degree seat opening facilitates easy transfers.
- Can be used over a toilet with adjustable seat-to-floor height, minimizing mess & maintaining client dignity.
- Locking casters improve transfer safety by preventing chair movement.
- Comfortable positioning armrests make clients feel safe & supported.
- Anti-slip footrest helps position clients securely in the chair.
- Front safety bar provides additional security during transport.
- Versatile use for transferring clients from bed, commode, showering area, or into other rooms.
- Adjustable seat height allows safe & easy transfers to various surfaces.
- 1 button electric lift design

FEATURES

- Soft to Touch Padded Seat, Arms, & Backrest
- Supportable Backrest
- Soft to Touch Coated Frame
- Waterproof Electric Push/Pull Lever
- Soft Padded Footplates
- Easy to Access Charging Box
- Safety Lock with Coated Plastic Cover
- 2 Locking Institutional Casters & 2 Non-locking Casters
- Padded Black Removable Cushion
- Seat Height Range: Approximately 16"–20"



REMOVABLE BLACK CUSHION