

BARIATRIC SOFT SEAT SHOWER COMMmode



700 SERIES

THE BARIATRIC SOFT SEAT SHOWER COMMmode CHAIR IS STRONG, DURABLE AND COMFORTABLE

The Bariatric Soft Seat Shower Commode with Footrest (700 Series) has a 700 pound weight capacity and a choice of 26", 28" or 30" seat width.

AVAILABLE IN 3 WIDTHS

FEATURES

- Soft seat included
- High density polyethylene flat seat
- Swing out arm rests
- User-friendly push-pull handle
- Backrest with upholstery and adjustable tension straps
- Slide out footrest & self storing
- Weight capacity 700 lbs
- 6 Deluxe Heavy Duty Casters: 5" x 1.25" (4 Brake and 2 Non-Brake)
- 10 qt commode



SPECS	Internal Width	External Width	Internal Depth	Seat Height	Overall Dimensions	Shipping Weight	Weight
BAR 726	26"	30"	23"	24"	Bar 726 45" H 30" W 28" D	36" x 36" x 51"	132lbs
BAR 728	28"	32"	23"	24"	Bar 728 45" H 32" W 28" D	36" x 36" x 51"	132lbs
BAR 730	30"	34"	23"	24"	Bar 730 45" H 34" W 28" D	36" x 36" x 51"	132lbs