


ELEVATE YOUR CONVENIENCE WITH A SUITCASE RAMP. NOW OFFERED IN LENGTHS OF 3", 4", 5" & 6" C/W TRAVEL BAG



EZ3RGT, EZ3, EZ4, EZ5 & EZ6

Ezee Life Single fold ramp Experience swift setup and convenience with our portable ramp featuring a single-fold design, perfect for scooters and wheelchairs. Constructed from lightweight aluminum, the SUITCASE Singlefold Ramp boasts an easy-to-carry design, complete with ergonomically-designed handles and a full-length hinge for added strength and stability.

Specifically engineered for high-traffic use, this durable ramp is offered in lengths ranging from 3', 4', 5' and 6', featuring an extruded slip-resistant tread for superior traction.



Key Features

- Extruded slip-resistant tread for superior traction.
- Wide enough for standard wheelchairs
- 600lbs (272kg) weight capacity
- Corrosion free marine grade aluminium material
- 2-way fold-away design for maximum convenience
- Carry handles and for easy transport
- Travel-Friendly Storage Bag for Ramps
- Raised safety rails to prevent roll-off
- Locking pins to secure your ramp when in use
- Certified to the highest international CE standards guaranteeing the quality and safety of each ramp (EN ISO 12100:2010)



SKU	PRODUCT NAME	OPEN SIZE (CM)	OPEN SIZE (INCH)	FOLD SIZE (CM)	FOLD SIZE (INCH)	PACKING SIZE (CM)	PACKING SIZE (INCH)	N.W (KG)	N.W (LBS)	BAG (KG)	BAG (LBS)	G.W (PRODUCT+BAG+CARTON) (KG)	G.W (PRODUCT+BAG+CARTON) (LBS)
EZ3RGT	3' Single Fold Portable Ramp with Grip Tape	72x92x5	28x36x2	36x92x8	14x36x3	40x94x10	16x37x4	7.7	17	0.9	2	9.6	21
EZ3R	3' Single Fold Portable Ramp	72x92x5	28x36x2	36x92x8	14x36x3	40x94x10	16x37x4	7.8	17	0.9	2	9.6	21
EZ4R	4' Single Fold Portable Ramp	72x122x5	28x48x2	36x122x8	14x48x3	40x125x10	16x49x4	10.2	22.5	1.1	2.4	12.3	27
EZ5R	5' Single Fold Portable Ramp	72x152.5x5	28x60x2	36x152.5x8	14x60x3	40x155x10	16x61x4	12.7	28	1.3	3	15.3	34
EZ6R	6' Single Fold Portable Ramp	72x183x5	28x72x2	36x183x8	14x72x3	40x185x10	16x73x4	15.2	34	1.5	3	18.2	40

EZ3RGT, EZ3, EZ5 & EZ6 Singlefold Ramp c/w with travel bag



34 Futurity Gate, Unit # 15 Concord Ontario L4K 1S6
 wk: 416-739-1267 | Toll Free: 1-888-634-5808 | fx 416-739-7171
 email: info@hpurehab.com | www.hpurehab.com