

EZee Life™

32" Folding Reacher

Model # CH3052



Limited Warranty:

Ezee Life™ products carry a one year warranty against manufacturing defects, faulty materials and workmanship only.

Any alterations, misuse, abuse, or accidental damage voids this warranty. Ezee Life™ reserves the right to replace or repair any that has become defective.

The warranty does not extend to consequential costs resulting from the fault or defect of the product: freight and travel costs, loss of earnings, or other expenses that one may incur.

Warranty is valid for the original purchaser only and the original receipt must be presented before any warranty options are considered. Please also be prepared to provide the date of purchase and serial number.

EZee Life™ Products

34 Futurity Gate, Unit # 15 Concord Ontario L4K 1S6



Description:

The EZee Life™ CH3052 reacher features a rotating head, magnetic tip molded jaw and trigger assembly. This reacher is a great item to help people who have bending or range of motion problems pick things off the floor or off high shelves safely.

Specifications:

- Jaw opening: 4"
- Length 32"
- Frame material: Aluminum tube
- Weight capacity: 18 lbs
- Product weight: 2 lbs

Features:

- 1/2" Aluminum shaft (ABS Handle & Jaws)
- 32" Length
- Folds for easy Storage
- Rotating head
- Magnetic tip
- Molded jaw and trigger for comfort
- Pressure sensitive jaw
- 4" Jaw opening

Instructions for Use:

To fold, grasp the hand grip and jaws and pull apart. Once ends are separated fold ends together.

To unfold, straighten out shaft, align ends together and insert jaw end into handle end.

Swivel jaws for more flexibility when picking up items.

Magnet located on the jaw tip to make picking up small metal objects easier.