

EZee Life™

Bed Rest Back Support

Model: CH3055



EZee Life™ Products

34 Futurity Gate, Unit 15
Concord, Ontario L4K 1S6



Limited Warranty:

Ezee Life™ products carry a one year warranty against manufacturing defects, faulty materials and workmanship only.

Any alterations, misuse, abuse, or accidental damage voids this warranty. Ezee Life™ reserves the right to replace or repair any part that has become defective.

The warranty does not extend to consequential costs resulting from the fault or defect of the product: freight and travel costs, loss of earnings, or other expenses that one may incur.

Warranty is valid for the original purchaser only and the original receipt must be presented before any warranty options are considered. Please also be prepared to provide the date of purchase and serial number.

Description:

The Ezee Life™ bed rest back support is a low cost alternative to a home hospital bed for those who don't require the features a regular hospital bed provides. Five back angle settings available.

Set up:

The Ezee Life™ bed rest back support comes fully assembled and ready to use. Simply set the back support to the desired back angle and place it on the bed.

Note – This product works best when it is placed against the headboard or wall to prevent slipping.

Specifications:

Back Width: 24"

Back Height: 25"

Frame Depth: 18"

Back Angles: 45°, 55°, 65°, 75° and 85°

Care and Cleaning:

To maintain the beauty of the original finish, clean with mild soap and rinse with warm water.

Wipe back support dry after use if required.

Precautions:

Periodically check assembly for continued tightness.

NEVER USE BACK SUPPORT IF ANY LOOSENESS IN THE ASSEMBLY IS NOTICED AS PERSONAL INJURY MAY RESULT.

Imported / Distributed by:

Factory Direct Medical
34 Futurity Gate, Unit # 15
Concord Ontario L4K 1S6
Local: 416-739-8393
Toll Free: 855-235-2400

